



# DOMESTIC PARTS LISTING

## Medical Equipment Washing & Decontamination Systems

Item	Description	UM	Status	Notes - Where Used	
<b>CHEMICAL INJECTOR PUMPS</b>					
2	225-22028	CHEM PUMP HOSE - YELLOW	Each	Std	Hose inside chemical pumps
2	225-36718	CHEM PUMP & MOTOR-BETA-115V,50/60HZ	Each	Std	Chem pump - Beta brand (green covers).
2	225-3618M	CHEM PUMP MOTOR ONLY	Each	Std	Motor only for Beta pump head
2	225-38850	CHEM PUMP REAR COVER: HSG ASSY	Each	Std	Chem pump - Beta brand (green covers).
2	225-51189	CHEM PUMP ROLLER ASSEMBLY	Each	Std	Chem pump - Beta brand (green covers).
2	225-59967	CHEM PUMP FRONT COVER: HSG Ø 0	Each	Std	Chem pump - Beta brand (green covers).
1	XM3-68-015	DETERGENT SUCTION HOSE-RED	Each	Std	
1	XM3-68-005	RINSE SUCTION HOSE-BLUE	Each	Std	
<b>CONTROL PANEL ASSEMBLY COMPONENTS (2004 forward)</b>					
1	256-21MD000BA5	SIEMENS LOGO CONTROLLER	Each	Std	Version 8
1	256-51HB000BA0	SIEMENS LOGO EXPANSION MODULE	Each	Std	Version 8
1	278-PS5R-C24	POWER SUPPLY 30 WATT 24VDC	Each	Std	Electrical Components
1	277-CLO0A310TJ	MOTOR CONTACTOR GE 3HP/230	Each	Std	Electrical Components
1	277-P9MEMHRN	SWITCH GE STOP BUTTON 40MM RED	Each	Std	Red Stop Button on Control Panel
1	277-P9MPLVGD	SWITCH GE START BUTTON - GREEN	Each	Std	Green Start Button
2	277-P9MSMD0N	SWITCH GE 2 POS BLK KNOB	Each	Std & Option	Run/Drain Switch and Lighted Wash Chamber Switch
2	277-P9MSMZ0N	SWITCH GE 3 POS BLK KNOB	Each	Std & Option	Wash Time & Wash Mode Switch
1 to 3	277-3M533	CAPCITOR MOUNTING KIT	Each	Std & Option	Electrical Components
1 to 3	277-97-F5705BX	CAPCITOR, RUN 5MFD 370 V	Each	Std & Option	Electrical Components
1	277-P9B01VN	CONTACT BLOCK GE 1 NC-STOP BTN	Each	Std	Electrical Components
1	277-P9B10VN	CONTACT BLOCK GE 1 NO - 3POS	Each	Std & Option	Electrical Components
8	238-9080GM6	TERM BLK W/C - SMALL BLOCK	Each	Std	Electrical Components
4	238-9080GR6	TERM BLK W/C - LARGE BLOCK	Each	Std	Electrical Components
<b>CONTROLS- SWITCHES, SENSORS, LEVEL CONTROLS</b>					
1	136-41535K91	PRESSURE SWITCH 1/4NPT 5-65PSI	Each	Option	Low Pressure Cntrl Switch (Get Nipple too)
1 or 2	235-2XC13	FLOAT SWITCH - HORIZONTAL LEVEL	Each	Std	Water Level High Float - Plastic Float Style
1	217-4C445	AIR BLOWER-TIMER DIGITAL AUTO OFF	Each	Option	
1	235-B11-204	LIMIT SWITCH-FR/DR 80-348 XTRN	Each	Std	Front door -white switch (lever type)
1	235-GSW-13	TOGGLE SWITCH 3POS CHEM PUMP	Each	Std	Chemical Pump Primer Switch
<b>CONTROL-FUSES</b>					
3 to 5	237-GMD 1.5A	FUSE 1.5 AMP GLASS TUBE	Each	Std	SA Drv Mtr, C-Pumps Inputs, Switches, Floats
1 to 3	237-6F057	FUSE 2 AMP GLASS TUBE	Each	Std	Chemical Pump Fuse
1	237-GMA 2.5	FUSE 2.5 AMP GLASS TUBE	Each	Std	PLC Outputs
1	237-GDC 5	FUSE 5 AMP GLASS TUBE	Each	Std or Option	Conveyor Motor
1	237-FNQ-R-5	FUSE 5 AMP CC MINI CRTRDGE (kldr-5)	Each	Std or Option	Transf (240v) Elect Sys (110v)-Drain pump
1	237-FRQ-R-8	FUSE 8 AMP CC MINI CRTRDGE (kldr-8)	Each	Std	Drain pump
1	237-FRN-R-2	FUSE 2 AMP FRN CARTTDGE	Each	Std	Control circuit-Drain Pump
1 or 2	237-FRN-R-20	FUSE 20 AMP FRN CRTRDGE (RFK)	Each	Std	Pump Motors & Main System Fuse
<b>DISINFECTANT APPLICATION SYSTEM-Standard and Optional</b>					
1	133-N18013208	PRESSURE GAUGE 0-60R 1/4NPT 2	Each	Option	Disinfectant Line Pressure Gauge
1	133-R43301N1LA	REGULATOR 3/8 NORGREN	Each	Option	Disinfectant Line Volume Regulator - 3/8"
1	138-202CT	DEMA VALVE INJECTOR 1/4" THREAD	Each	Option	Disinfectant Chemical Injector - 1/4"
4 to 8'	138-100 12	DEMA VINYL TUBING	Lin Ft.	Option	Disinfectant Chemical Tubing
1	138-24 11MP	DEMA FOOT STRAINER	Each	Option	Disinfectant Chemical Tubing Hose Strainer
1	151-50775K175	ADP POL BULKHEAD 3/8 POLY	Each	Option	Disinfect line through chamber wall
1	151-5532K644	ADP POL EL 3/8POLY X 3/8MNPT	Each	Option	Disinfect line to disinfect SA header
4	151-5532K422	ADP POL STRT 1/4POLY X 1/4MNPT	Each	Option	Disinfect header & rod style level probe
1	151-5532K944	ADP POL STRT 3/8POLY X 3/8FNPT	Each	Option	Disinfect Dema valve out line
4	153-C9PMN-BB	ADP CHR M EL 1/4POLY X 1/4MNPT	Each	Option	Disinfectant Chemical Line to Bottom Corner SA
<b>FILTERS &amp; SCREENS</b>					
1	M3-30-300	RESERVOIR FILTER SCREEN -3K SQUARE	Each	Std	
1	170-10373	FILTER SCREEN, FILL PORT 1	Each	Std	Filter Screen located inside Fill Line Hose (washer style)
1	AQ-170	WATER FILTRATION ASSM 1/2 NPT	Assembly	Option	Housing W/Plumbing & 2 10 mic carts
1	141-15839-75	FILTER CART HOUSING 1/2" NPT (CART 10 MIC FILTS)	Each	Option	Without Plumbing Fittings



# DOMESTIC PARTS LISTING

## Medical Equipment Washing & Decontamination Systems

Item	Description	UM	Status	Notes - Where Used	
1	141-EPM-10	FILTER CART-10 MICRONS	Each	Option	0
<b>GASKETS</b>					
1	M3-75-001	FRONT DOOR GASKET KIT 3K	Each	Std	143
1 or 2	XM4321-75-085	FRONT DOOR WINDOW GASKET	Each	Option	Front Window Gasket
<b>GFCI</b>					
1	234-54890R	GFCI-120V, 20 AMP NO POWER CORD	Each	Std	
1	234-2000-0050	GFCI-120V, 20 AMP WITH POWER CORD	Assembly	Std	
<b>MISC</b>					
1	291-LAM 100ARS	LITE BULB -FOR CHAMBER FIXTURE	Each	Option	Rough Service Bulb -Lighted Chamber Opt.
1	340-6X10-1/8	GLASS- LIGHTED CHAMBER WINDOW	Each	Option	Lighted Wash Chamber Replacement Window
1	340-12X12X.030	GLASS-FRONT DOOR WINDOW	Each	Option	Safety Glass Replacement Window
4	321-6466K27	FRONT WHEEL CARRIAGE BEARING-PVC	Each	Option	Select PVC or SS Bearing for replacement
4	321-590BK18	FRONT WHEEL CARRIAGE BEARING-SS	Each	Option	Select PVC or SS Bearing for replacement
1	170-10031640	FILL OR DRAIN HOSE 10-FOOT LONG	Each	Std	Vulcanize hose ends 180 degree rating.
1	AQ-195	AQ OPERATIONS MANUAL	Each	Std	AQ-4000 Operations Manual - Printed & Binded
<b>NOZZLES</b>					
23	161-NBC1565	NOZZLE 1/4 1.5 GPM 60D FULL CONE	Each	Std	Stairionay Arm Units Only.
<b>PUMPS &amp; MOTORS</b>					
1	228-EBACDU200/115D1C	PUMP & MOTOR ASSEMBLY 1.5 HP EBARA	Assembly	Option	
		PUMP & MOTOR ASSEMBLY 1 HP EBARA	Assembly	Std	
1	211-MAGB719	PUMP MOTOR ONLY 1HP 120V 60HZ 56J-FRAME	Each	Std	Domestic 1K, 2K, & 3K
		PUMP MOTOR ONLY 1.5 HP 120V 60HZ 56J-FRAME	Each	Option	Upgrades, 3K, 3.7K & 8k Models
1	EB6625001	PUMP MECH SEAL	Assembly	Std	Internal Pump Mechanical Seal
1	222-EB165X3	PUMP O-RING EBARA CDU/70/5NB-2H	Each	Std	ALL PUMPS
1	EBCDU120/1NB	PUMP - PUMP HEAD ONLY- KIT, 1 HP 120/1	Each	Std	Domestic 250 1K 1.5K, 2K, 2.7K
1	EBCDU200/1NB	PUMP - PUMP HEAD ONLY- KIT, 1.5 HP 120/1	Each	Option	Domestic 3K, 3.7K, 8k
1	222-251450788	PUMP IMPELLER ONLY 1 HP EBARA CDU 120/1	Each	Std	Domestic 1K, 2K, & 3K
1	0	PUMP IMPELLER ONLY 1.5 HP EBARA CDU 120/1	Each	Option	
1	222-340280085	PUMP MOUNTING FEET	Each	Std	All pumps
1	228-3201-96	DRAIN- PUMP & MOTOR-AMT-ALUM. CAST	Assembly	Std	Cast Alum. Pump head 1-1/4 x 1-1/4
1	217-4C445	AIR BLOWER MOTOR 495 CFM	Each	Option	Small AQ Units up to 37K
1 TO 3	211-4Z062	GEAR MOTOR 14 RPM AC - 4K SA	Each	Option	Spray Arm Drive Motor
<b>SPRAY ARMS</b>					
8 TO 12	317-9452K39	O-RING BUNA DASH 219 B/HD 1/8	Each	Std or Option	Interior bulkhead O-Ring for spray arms
1 TO 4	M3-35-035	SPRAY ARM BEARING SUPPORT BRKT	Assembly	Option	FOR OLDER AQ UNITS
4 TO 8	M3-35-005	SA BEARING L/R SIDE	Each	Std	UHMW Material - SA Bearing - both Left & Right Sides
<b>THERMOMETERS</b>					
1	316-3946K125	THERMOMETER RND DIAL 3" 1/2NPT	Each	Option	Fill Line Thermometer - Optional
1	316-38675K13	THERMOMETER *C* CONTROL PANEL	Each	Option	International, option
1	316-4141K11	THERMOMETER ECONOMY PANEL	Each	Option	Control Panel Water Temperature Display - Optional
<b>TRACK</b>					
2	M32-55	BRUSH ASSEMBLY FRONT 4" - EACH	Assembly-Each	Option	Track, requires two
2	M32-56	BRUSH ASSEMBLY REAR 6" - EACH	Assembly-Each	Option	Track, requires two
3	M3-49-065	TRACK DRIVE BAR SHAFT-3K	Each	Option	Drive bars for wheel rotation (3)
4	321-1F636	PILLOW BLOCK BEARING 5/8 W SS HOUSING	Each	Option	
3	M3-49-100	TRACK SPROCKET-SMALL 14 1/2	Each	Option	Small sprockets on ends of drive bars (3)
2	M3-49-110	TRACK SPROCKET-LARGE 40 5/8	Each	Option	Large sprockets on ends of front drive bar (2)
2	5M125-0	TRACK COTTER PIN-3/16 X 3 STAINLESS	Each	Option	Large sprocket pins
3	98355A160	TRACK COTTER PIN-1/8 X2-1/4 STAINLESS	Each	Option	Small sprocket pins
1	M3-49-200	TRACK CHAIN-LEFT PLAST 31-5/8	Assembly	Option	Plastic chain on left side
1	M3-49-210	TRACK CHAIN-RIGHT PLAS 25-5/8	Assembly	Option	Plastic chain on right side
1	M3-49-220	TRACK CHAIN-MOTOR SS 54-11/16	Assembly	Option	Plastic/SS chain from drive motor to large sprock
1	M3-159	CHAIN OPENING SEAL	Each	Option	
1	324-35PC CL	CHAIN CONNECTING LINK	Each	Option	Track drive chain - SS links
0	324-FGSATL3500	CHAIN-DELRLN-35 PITCH	Lin. Ft.	Option	



# DOMESTIC PARTS LISTING

## Medical Equipment Washing & Decontamination Systems

Item	Description	UM	Status	Notes - Where Used	
0	324-FT35PC	CHAIN-SS LINK-35 PITCH	Lin. Ft.	Option	
4	321-6389K552	BEARING SLEEVE, FLNGD, NYLON, 5/8	EA	Option	Track drive bar (2) & door handle (2)
2	321-11951-SS	DRIVE ROD END JOINT-SS	Each	Option	
6	322-6436K136	COLLAR-ALUM SET 5/8 BORE 2PC	Each	Option	Track drive bars set collars (4)
1	323-35BF101/2	DRIVE MOTOR SPROCKET 10 TOOTH FORGED 1/2"	Each	Option	Track drive motor sprocket (1)
TUBING, HOSES, FITTINGS (AND ACCESSORIES)					
0	151-5532K422	ADP POL STRT 1/4POLY X 1/4MNPT	Each	Option	Adaptor for Data Collections System
4	153-C9PMN-BB	ADP CHROME EL 1/4 MNPT X 1/4 POLY	Each	Option	Disinfect Nozzels of Bottom Corner Spray Arms
0	151-5532K444	ADP POL STRT 3/8POLY X 3/8MNPT	Each	Option	
0	151-5532K488	ADP POL STRT 3/8POLY X 1/4MNPT	Each	Option	
1	151-5532K644	ADP POL EL 3/8POLY X 3/8MNPT	Each	Option	
1	151-5532K822	ADP POL TEE 1/4POLY X 1/4MNPT	Each	Option	Chem hose connect at reservoir
2'	172-0100	POLY TUBE 1/4 RED	Lin. Ft.	Std	Chemical Detergent Suction Hose
2'	172-0120	POLY TUBE 1/4 BLUE	Lin. Ft.	Std	Chemical Rinse Suction Hose
40'	172-0140	POLY TUBE 1/4 NATURAL	Lin. Ft.	Std	Chemical Injection Hose - Into Wash Chamber
2'	172-0340	POLY TUBE 3/8 NATURAL	Lin. Ft.	Option	In Wash Chamber & from Dema device.
2	170-10036040	FILL OR DRAIN HOSE 10'	Each	Std	
Varies	170-HPD-6 41662	HOSE 1 BLACK	Lin. Ft.	Std	Spray arm distrib., drain pump suction
4'	170-HPD-9 61662	HOSE 1-1/2 BLACK	Lin. Ft.	Std	Main pump suction
0	170-102-1000	HOSE 3/4 BLACK	Lin. Ft.	Not-Std	Drain pump suction
2	170-70815T43	HOSE CONNECT 1/2 GARDEN FEMALE	Each	Std	Fill and Drain connections