wheelchairs and other hard-to-clean medical equipment.

# AQ-F250 Features & Options

## COMPARE PRICING WITH OUR MOST **BASIC WASHER--**THE AQ-F250 SPRITE

#### **Standard Features:**

- Double Wall Stainless Steel Chamber
- PVC Spray Arms
- Variable Wash Time Cycles (3, 5, or 7 Minute)
- Wash and Rinse Cycles
- Automatic Chemical Injection
- LCD Status Display
- Heavy-Duty Caster Wheels
- Fill & Drain Hose Set
- 1 HP Motor/ Stainless Steel Pump
- 115 Volt 20 Amp AC Operation

### **Optional Features:**

- Stainless Steel Spray Arm Upgrade
- Front Door Window



#### **Total Outer Dimensions**

• Width 33-1/2"

• Length 64"

• Height 62"

**Effective Chamber Dimensions** 

• Width 26-1/2"

• Length 51-1/4"

• Height 41-1/2"

Like all Aqua Phase models, in order to operate the machine once the unit is located, all you need to do is hook up two hoses — one for the hot water fill and one for the drain, connect the Aqua Phase to an appropriate power source and begin saving time, improving efficiency and enhancing your infection control protocol.