



Vashers and Sanitizers to Automatically Wash and Apply Disinfectant to Wheelchairs and Other Hard-to-Clean Medical Equipment.

# **TAKE-DOWN MODEL**

# AQ-4200 model is slightly narrower than the AQ-4000

## STANDARD FEATURES

### Structural

- Double Wall Stainless Steel Chamber
- Stainless Steel Structural Frame
- Stainless Steel Control Panel Enclosure
- Oscillating Stainless Steel Spray Arms
- Chamber Floor Deck & Loading ramps
- Caster Wheels & Leveling Feet
- Seal Welded Wash Chamber
- Front Door 12ga Stainless Steel
- Front Door Safety Window
- Disinfectant Application System
- Lighted Wash Chamber
- Fill Line Temperature Gauge 3" Dial

### Controls & Electrical

- Programmable Logic Control Syst (PLC)
  Capacity Gallons/min: 90
- 24 Volt DC Control System
- LCD Status Display
- Level Control Self Diagnostics
- Wash Cycle Pressure Monitoring
- Non-Reset Cycle Counter
- GFCI Built in Circuit
- 208/230 Volt, 25 Amp, 50/60 Hertz SP
- UL, CSA, CE/EU

# **O**PTIONAL FEATURES & **U**PGRADES -

- Utility Wash Tray Set Stainless Steel
- Digital Temperature Display
- Pass-Through Design (Entry/Exit Doors)
- Conveyor Belt Transfer System
- Data Collections Computer System
- Fire Equipment/KEDS Wash Rack Stainless Steel
- Chamber Swing Arms to wash Backboards

### Exterior

- Length: 98-7/16"
- Width: 63-1/2"
- Height: 73"

- **DIMENSIONS & WEIGHT** —

- Chamber with Floor Deck • Length: 93"
  - Width: 38-3/4"
  - Height: 47-5/16"

- Automatic Filling & Draining Mechanical
  - 3 HP Motor/Stainless Steel Pump

### Wash System Nozzle Data

Automatic Chemical Injection

- Wash & Rinse System Nozzles: 60
- Disinfectant Nozzles: 36
- Wash PSI: 28
- Capacity Gallons/cycle: 630
- Capacity Gallons/hour: 5400

### Extras

- Fill & Drain Hose Set 10' Each
- Loading Ramp Set (2 @ 48"L x 16"W)
- Indicate Door Swing(s) to fit your location

  - Disinfectant Application Spray Wand (Exterior Mounted)
  - Wash-Down Spray Gun & Hose (Exterior Mounted)
  - Back Flow Prevention Devices / System
  - Water Supply Temperature Mixing Valve (1/2" or 3/4" Line Size)

AQ-4200 Model shown with a painted Control Panel Enclosure

- Disinfectant Mixing Station (Exterior Mounted)
- Specialized Cycle Sequencing
- Air Dry Blower
- Crate/Pallet Dimensions
  - Length: 112"
  - Width: 84"
  - Height: 84"

 2175 lbs. • 986 Kilo

Weight

Approximate Shipping

V130426

Door Opening Requirements: Width: 47-3/4" x Height: 67-3/4"

Once the Aqua Phase unit is in its location, all you need to do is hook up two hoses - one for the hot water fill line and one for the drain line. Then plug the Aqua Phase in an appropriate power source & begin saving time, improving efficiency & enhancing your infection control protocol.

> Aqua Phase, A Division of Mid-State Stainless, Inc. / 330 W. Benson Ave. / PO Box 228, Grantsburg, WI 54840 / USA p: 800-208-9274 | f: 715-463-5853 | aquaphase@medwashers.com | www.aqua-phase.com





- Clean Only, Clean & Apply Disinfectant,
  - and Apply Disinfectant Only modes